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Wang, J.; Chow, W.; Leung, D.; Chang, J. Application of ultrahigh-performance liquid chromatography and electrospray ionization quadrupole orbitrap high-resolution mass spectrometry for determination of 166 pesticides in fruits and vegetables. J. Agric. Food Chem. 2012, Dec 12;60(49):12088-104.

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Wang, J.; Chow, W.; Chang, J.; Wong, J. W. Ultrahigh-performance liquid chromatography electrospray ionization Q-Orbitrap mass spectrometry for the analysis of 451 pesticide residues in fruits and vegetables: method development and validation. J. Agric. Food Chem. 2014, Oct 22;62(42):10375-91.

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Zomer, P. and Mol H. G. J. Simultaneous quantitative determination, identification and qualitative screening of pesticides in fruits and vegetables using LC-Q-Orbitrap™-MS. Food Addit. Contam: Part A. 2015 Sep 15:1-9. [Epub ahead of print].

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Kaufmann, A.; Butcher, P.; Maden, K.; Walker, S.; Widmer, M. Multi-residue quantification of veterinary drugs in milk with a novel extraction and clean up technique: salting out supported liquid extraction (SOSLE). Anal. Chim. Acta. 2014, Apr 11;820:56-68.

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Jia, W.; Chu, X.; Ling, Y.; Huang, J.; Chang, J. Multi-mycotoxin analysis in dairy products by liquid chromatography coupled to quadrupole orbitrap mass spectrometry. J. Chromatogr. A. 2014, Jun 6;1345:107-14.

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Dzuman, Z.; Zachariasova, M.; Veprikova, Z.; Godula, M.; Hajslova J. Multi-analyte high performance liquid chromatography coupled to high resolution tandem mass spectrometry method for control of pesticide residues, mycotoxins, and pyrrolizidine alkaloids.

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Jia, W.; Chu, X.; Ling, Y.; Huang, J.; Chang, J. High-throughput screening of pesticide and veterinary drug residues in baby food by liquid chromatography coupled to quadrupole Orbitrap mass spectrometry. J. Chromatogr. A. 2014, Jun 20;1347:122-8.

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Gomez-Perez, M. L.; Romero-Gonzalez, R.; Martinez Vidal, J. L.; Garrido Frenich, A. Analysis of veterinary drug and pesticide residues in animal feed by high-resolution mass spectrometry: comparison between time-of-flight and Orbitrap.

Food Addit. Contam. Part A Chem. Anal. Control Expo. Risk. Assess. 2015, Mar 18:1-10.

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Peer reviewed journals:

LC and LC-MS methods

Title	Authors	Publication
Release of non-intentionally added substances (NIAS) from food contact polycarbonate: Effect of aging	Chiara Bignardi, Antonella Cavazza, Carmen Laganà, Paola Salvadeo, Claudio Corradini	Food Control, Volume 71, January 2017, Pages 329–335
UHPLC-high-resolution mass spectrometry determination of bisphenol A and plastic additives released by polycarbonate tableware: influence of aging and surface damage	Chiara Bignardi, Antonella Cavazza, Carmen Laganà, Paola Salvadeo, Claudio Corradini	Analytical and Bioanalytical Chemistry, October 2015, Volume 407, Issue 26, Pages 7917–7924
Targeted and untargeted data-dependent experiments for characterization of polycarbonate food-contact plastics by ultra-high-performance chromatography coupled to quadrupole Orbitrap tandem mass spectrometry	Chiara Bignardi, Antonella Cavazza, Claudio Corradini, Paola Salvadeo	Journal of Chromatography A, Volume 1372, 12 December 2014, Pages 133–144
Liquid chromatography–full scan-high resolution mass spectrometry-based method towards the comprehensive analysis of migration of primary aromatic amines from food packaging	Monica Mattarozzi, Francesca Lambertini, Michele Suman, Maria Careri	Journal of Chromatography A, Volume 1320, 13 December 2013, Pages 96–102

Methods

Title	Authors	Publication
The challenge of identifying non-intentionally added substances from food packaging materials: A review	C. Nerin, P. Alfaro, M. Aznar, C. Domeño	Analytica Chimica Acta, Volume 775, 2 May 2013, Pages 14–24
Rapid qualitative analysis of phthalates added to food and nutraceutical products by direct analysis in real time/Orbitrap mass spectrometry	Randy L. Self, Wen-Hsin Wu	Food Control, Volume 25, Issue 1, May 2012, Pages 13–16
PCDD and PCDF levels in paper with food contact	H. Beck, A. Droß, W. Mathar	Chemosphere, Volume 25, Issues 7–10, October–November 1992, Pages 1533–1538

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Characterizing unknowns in food packaging using GC Orbitrap Mass Spectrometry

Dominic Roberts¹, Jesus Varela², Yves-Alexis Hammel² and Paul Silcock¹ ¹ Thermo Fisher Scientific, Runcorn, UK ² Nestle, Lausanne, Switzerland, Thermo Scientific application note: AN 10527

352 RAJSKI ET AL.: JOURNAL OF AOAC INTERNATIONAL VOL. 101, NO. 2, 2018

SPECIAL GUEST EDITOR SECTION

Coupling Ion Chromatography to Q-Orbitrap for the Fast and Robust Analysis of Anionic Pesticides in Fruits and Vegetables

LUKASZ RAJSKI, FRANCISCO JOSÉ DÍAZ GALIANO, VÍCTOR CUTILLAS, and AMADEO R. FERNÁNDEZ-ALBA¹



University of Almería, Department of Chemistry and Physics, Agrifood Campus of International Excellence (CeiA3), European Union Reference Laboratory for Pesticide Residues in Fruits and Vegetables, 04120 Almería, Spain

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thermoscientific

APPLICATION NOTE

Characterizing unknowns in food packaging using GC Orbitrap Mass Spectrometry

Authors: Dominic Roberts¹, Jesus Varela², Yves-Alexis Hamme³ and Paul Silcock¹
¹Thermo Fisher Scientific, Runcorn, UK
²Nestlé, Lausanne, Switzerland

Key Words

Food packaging, Q Exactive GC, Orbitrap mass spectrometry, unknown identification, structural elucidation, food safety

Analytica Chimica Acta 760 (2013) 65–73

Targeted analysis with benchtop quadrupole-orbitrap hybrid mass spectrometer: Application to determination of synthetic hormones in animal urine

No. 10527

CrossMark
Click for updates

Food Additives & Contaminants: Part A

Publication details, including instructions for authors and subscription information: <http://www.tandfonline.com/loi/tafca20>

Analysis of veterinary drug and pesticide residues in animal feed by high-resolution mass spectrometry: comparison between time-of-flight and Orbitrap

María Luz Gómez-Pérez², Roberto Romero-González², José Luis Martínez Vidal¹ & Antonia Garrido Frenich¹

¹Research Group "Analytical Chemistry of Contaminants", Department of Chemistry and Physics, Research Centre for Agricultural and Food Biotechnology (BITAL), University of Almería, Agrifood Campus of International Excellence, Almería, Spain
 Published online: 18 Mar 2015.

JOURNAL OF AGRICULTURAL AND FOOD CHEMISTRY

Article
pubs.acs.org/JAF

Ultrahigh-Performance Liquid Chromatography Electro spray Ionization Q-Orbitrap Mass Spectrometry for the Analysis of 451 Pesticide Residues in Fruits and Vegetables: Method Development and Validation

Jian Wang^{a,b}, Willis Chow,^c James Chang^d and Jon W. Wong^e

^aCanadian Food Inspection Agency, Calgary Laboratory, 3650-36th Street N.W., Calgary, Alberta T2L 2L1, Canada

^bThermoFisher Scientific, 355 River Oaks Parkway, San Jose, California 95134, United States

^cCenter for Food Safety and Applied Nutrition, U.S. Food and Drug Administration, 5100 Paint Branch Parkway, College Park, Maryland 20740, United States

DOI: 10.1002/anie.201401579

PAPER IN FOREFRONT

Liquid chromatography Orbitrap mass spectrometry with simultaneous full scan and tandem MS/MS for highly selective pesticide residue analysis

María del Mar Gómez-Ramos¹, Lukasz Rajski¹, Horacio Heinzen², Amadeo R. Fernández-Alba¹

352 RAJSKI ET AL.; JOURNAL OF AOAC INTERNATIONAL VOL. 101, NO. 2, 2018

SPECIAL GUEST EDITOR SECTION

Coupling Ion Chromatography to Q-Orbitrap for the Fast and Robust Analysis of Anionic Pesticides in Fruits and Vegetables

LUKASZ RAJSKI, FRANCISCO JOSÉ DÍAZ GALIANO, VÍCTOR CUTILLAS, and AMADEO R. FERNÁNDEZ-ALBA¹
 University of Almería, Department of Chemistry and Physics, Agrifood Campus of International Excellence (CeiA3), European Union Reference Laboratory for Pesticide Residues in Fruits and Vegetables, 04120 Almería, Spain

Analytica Chimica Acta 820 (2014) 56–68

Contents lists available at ScienceDirect

Analytica Chimica Acta

journal homepage: www.elsevier.com/locate/aca

Multi-residue quantification of veterinary drugs in milk with a novel extraction and cleanup technique: Salting out supported liquid extraction (SOSLE)

A. Kaufmann^a, P. Butcher, K. Maden, S. Walker, M. Widmer

Official Food Control Authority of the Canton of Zurich, Fehrenstrasse 15, Zurich 8032, Switzerland

CrossMark

thermoscientific
WHITE PAPER 65147
Selecting the best Q Exactive Orbitrap mass spectrometer scan mode for your application

reviewed journals:

LC and LC-MS methods

Title	Authors	Publication
Release of non-intentionally added substances (NIAS) from food contact polycarbonate: Effect of aging	Chiara Bignardi, Antonella Cavazza, Carmen Laganà, Paola Salvedo, Claudio Corradini	Food Control, Volume 71, January 2017, Pages 329–335
LHPLC-high-resolution mass spectrometry determination of bisphenol A and plastic additives released by polycarbonate: influence of aging and surface damage	Chiara Bignardi, Antonella Cavazza, Carmen Laganà, Paola Salvedo, Claudio Corradini	Analytical and Bioanalytical Chemistry, October 2015, Volume 407, Issue 26, Pages 7917–7924
Targeted and untargeted data-dependent experiments for characterization of polycarbonate food-contact plastics by ultra-high performance chromatography coupled to quadrupole-Orbitrap tandem mass spectrometry	Chiara Bignardi, Antonella Cavazza, Claudio Corradini, Paola Salvedo	Journal of Chromatography A, Volume 1372, 12 December 2014, Pages 133–144
Liquid chromatography–full scan–high resolution mass spectrometry-based method towards the comprehensive analysis of migration of primary aromatic amines from food packaging	Monica Mattarozzi, Francesca Lambertini, Michele Sumari, Maria Careri	Journal of Chromatography A, Volume 1320, 13 December 2013, Pages 96–102

Methods

Title	Authors	Publication
The challenge of identifying non-intentionally added substances from food packaging materials: A review	C. Nerin, P. Alfaro, M. Aznar, C. Domeño	Analytica Chimica Acta, Volume 775, 2 May 2013, Pages 14–24
Rapid qualitative analysis of zwitterions added to food and nutraceutical products by direct analysis in real time/Orbitrap mass spectrometry	Randy L. Seif, Win-Hsin Wu	Food Control, Volume 25, Issue 1, May 2012, Pages 13–16
PCDD and PCDF levels in paper with food contact	H. Beck, A. Droll, W. Mether	Chemosphere, Volume 25, Issues 7–10, October–November 1992, Pages 1533–1538

Journal of Bioanal Chem

DOI: 10.1007/s00216-015-8709-z

PAPER IN FOREFRONT

Liquid chromatography Orbitrap mass spectrometry with simultaneous full scan and tandem MS/MS for highly selective pesticide residue analysis

María del Mar Gómez-Ramos¹, Lukasz Rajski¹, Horacio Heinzen², Amadeo R. Fernández-Alba¹

Journal of Chromatography A, 2014 (2014) 107–114
 Contents lists available at ScienceDirect
 Journal of Chromatography A
 journal homepage: www.elsevier.com/locate/chroma

Multi-mycotoxin analysis in dairy products by liquid chromatography coupled to quadrupole orbitrap mass spectrometry

Wei Jia^{a,b}, Xiaogang Chu^{a,b,c}, Yun Ling^b, Junrong Huang^d, James Chang^c

^a College of Chemistry & Chemical Engineering, Shaanxi University of Science & Technology, Xi'an 710021, China

^b Institute of Food Safety, Chinese Academy of Inspection and Quarantine, Beijing 100123, China

^c ThermoFisher Scientific, 355 River Oaks Parkway, San Jose, CA 95134, USA