
MALDI Biotyper

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Pál Róbert

Flextra-Lab Kft.

- **> \$1.7 billion revenue p.a. (2017)**
- **Founded in 1960 in Karlsruhe (DE)** by Prof. Günther Laukien with Nuclear Magnetic Resonance products

- Worldwide sales and service network with over 6.000 employees
- Global market leadership in many segments

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
MALDI Biotyper Timeline

- 1960 Foundation of the Bruker company
- 1988 Establishment of the mass spectrometry business in Bremen, Germany
- 1991 Introduction of the **first MALDI-TOF** mass spectrometer
- 2004 Launch of the MALDI Biotyper as a research tool
- **2008** **First MALDI Biotyper systems in routine laboratories**
- 2009 Launch of the MALDI Biotyper IVD system
- 2013 FDA Clearance of the MALDI Biotyper CA system
- **12/2017** **AOAC-OMA* approvals, OMA#2017.09 & OMA#2017.10**
- **4/2018** **> 2,900 MALDI Biotyper systems sold**
- **2/2018** **ISO 16140-part 6 validation by MicroVal, recognized by the EU 2073 as an alternative to confirm specific pathogens**

*AOAC-OMA: American Organisation of Analytical Chemists - Official Method of Analysis

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MALDI Biotyper Mass Spectrometer

- Benchtop instrument
- Unattended Operation

Dedicated MBT Compass Software for

- Automated data acquisition
- Automated data processing
- Automated pattern matching
- Report creation

MBT Compass Reference Library

- 2,748 species entries
- Continuously maintained & updated

MBT Subtyping Module



- Fast & easy *Listeria* species differentiation
- Resistance marker detection

Consumables

- Bruker HCCA Matrix
- Bruker Bacterial Test Standard
- MBT Biotarget 96


MALDI Biotyper - Basics

Matrix **A**ssisted **L**aser **D**esorption / **I**onization
Time **O**f **F**light Mass Spectrometry



MALDI Biotyper - Basics

Matrix **A**ssisted **L**aser **D**esorption / **I**onization
Time **O**f **F**light Mass Spectrometry


MALDI Biotyper - Basics

Matrix **A**ssisted **L**aser **D**esorption / **I**onization
Time **O**f **F**light Mass Spectrometry

MALDI Biotyper - Basics

Robust identification method,
as it relies on highly abundant proteins


E. coli

**Mass range calibrated:
2,000 – 20,000 Da**

MALDI Biotyper

Workflow independent from culture medium


Pseudomonas oleovorans

Species ID does not depend on growth conditions, or medium selection, as tested during the certification studies presented later

MALDI Biotyper

1 System – 1 Workflow: summary

Short time-to-result,
including

- Sample transfer
- Bacterial Test Standard for QC
- Drying of the HCCA matrix
- Mass spectrum acquisition
- Spectra matching

N° of isolates	TTR
1 isolate	~ 15 min
48 isolates	~ 35 min
95 isolates	~ 60 min


Same workflow for bacteria, yeasts and molds

5 to 10% alternative sample prep methods,
i.e. Extended Direct Transfer or Extraction method

MALDI Biotyper - Bioinformatics

Unsupervised Main Spectra (MSP) concept

Score base pattern matching

Unknown microorganism is matched against each Main Spectrum in the library

Calculation of a matching score based on:

Matches MSP to unknown

% matches of the reference spectrum (e.g. 6/10)

Matches unknown to MSP


% matches of the unknown spectrum (e.g. 6/20 = 3/10)

Correlation of intensities

value of intensity correlation

Range	Interpretation	Symbols	Color
2.000 - 3.000	High Confidence Identification	(+++)	green
1.700 - 1.999	Low Confidence Identification	(+)	yellow
0.000 - 1.699	No Organism Identification Possible	(-)	red

MALDI Biotyper RUO library Updates since 2010

Easy validation of library updates

MBT Filamentous Fungi Library

152 species / species groups

152 species / species groups

<i>Absidia caerulea</i>	<i>Aspergillus tamarit</i>	<i>Exophiala dermatitidis</i>
<i>Absidia glauca</i>	<i>Aspergillus terreus</i>	<i>Fusarium aquaeductum</i>
<i>Acaulium acremonium</i>	<i>Aspergillus tritici</i>	<i>Fusarium avenaceum</i>
<i>Alternaria alternata</i>	<i>Aspergillus unguis</i>	<i>Fusarium cerealis_culmorum_group</i>
<i>Arthrinium arundinis</i>	<i>Aspergillus ustus</i>	<i>Fusarium chlamydosporum</i>
<i>Arthrinium phaeospermum</i>	<i>Aspergillus versicolor</i>	<i>Fusarium delphinoides</i>
<i>Arthroderma amazonicum</i>	<i>Aspergillus westerdijkiae</i>	<i>Fusarium dimerum</i>
<i>Arthroderma eboreum</i>	<i>Aureobasidium melanogenum</i>	<i>Fusarium equiseti</i>
<i>Arthroderma flavescens</i>	<i>Aureobasidium pullulans</i>	<i>Fusarium incarnatum</i>
<i>Arthroderma gloriae</i>	<i>Beauveria bassiana</i>	<i>Fusarium oxysporum</i>
<i>Arthroderma lenticulare</i>	<i>Botrytis cinerea</i>	<i>Fusarium petrophilum</i>
<i>Arthrographis kalrae</i>	<i>Byssochlamys fulva</i>	<i>Fusarium proliferatum</i>
<i>Aspergillus calidoustus</i>	<i>Byssochlamys nivea</i>	<i>Fusarium solani</i>
<i>Aspergillus clavatus</i>	<i>Byssochlamys spectabilis</i>	<i>Fusarium sp</i>
<i>Aspergillus flavus_oryzae_group</i>	<i>Chaetomium globosum</i>	<i>Fusarium verticillioides</i>
<i>Aspergillus fumigatus</i>	<i>Chrysosporium keratinophilum</i>	<i>Isaria farinosa</i>
<i>Aspergillus glaucus</i>	<i>Chrysosporium shanxiense</i>	<i>Lasiodiplodia sp</i>
<i>Aspergillus izukae</i>	<i>Cladosporium cladosporioides</i>	<i>Lichtheimia corymbifera</i>
<i>Aspergillus japonicus</i>	<i>Cladosporium herbarum</i>	<i>Lomentospora prolificans</i>
<i>Aspergillus lentulus</i>	<i>Cladosporium sp</i>	<i>Metarhizium marquandii</i>
<i>Aspergillus minisclerotigenes</i>	<i>Cladosporium sphaerospermum</i>	<i>Microascus melanosporus</i>
<i>Aspergillus montevidensis</i>	<i>Coniochaeta hoffmannii</i>	<i>Microsporium canis</i>
<i>Aspergillus nidulans</i>	<i>Coniochaeta mutabilis</i>	<i>Monascus ruber</i>
<i>Aspergillus niger</i>	<i>Cunninghamella elegans</i>	<i>Monilinia laxa</i>
<i>Aspergillus nomius</i>	<i>Curvularia clavata</i>	<i>Mucor circinelloides</i>
<i>Aspergillus ochraceus</i>	<i>Curvularia lunata</i>	<i>Nannizzia fulva</i>
<i>Aspergillus parasticus</i>	<i>Curvularia pallenscens</i>	<i>Nannizzia gypsea</i>
<i>Aspergillus penicillioides</i>	<i>Curvularia verruculosa</i>	<i>Nannizzia incurvata</i>
<i>Aspergillus pseudoglaucus</i>	<i>Dichotomopilus funicola</i>	<i>Nannizzia persicolor</i>
<i>Aspergillus pulvinus</i>	<i>Didymella aurea</i>	<i>Nannizzia praecox</i>
<i>Aspergillus sclerotiorum</i>	<i>Epicoccum nigrum</i>	<i>Paraphyton cooki</i>
<i>Aspergillus sp[4]</i>	<i>Epicoccum sorghinum</i>	<i>Penicillium brevicompactum</i>
<i>Aspergillus sydowii</i>	<i>Epidemophyton floccosum</i>	<i>Penicillium chrysogenum</i>

152 species / species groups, for which contributions have been received from over 20 laboratories across 8 countries

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
MALDI Biotyper

Connect to online available libraries

e.g. MicrobeNet from CDC Atlanta

- Free available
- Internet connection in the lab needed
- Possible because of the MBT open microbiology concept

- Status May 2018:

1004 MSPs, all on rare and unusual species, to complement the Bruker MBT Compass Library, continuous expansion

- Future Expansions: rare pathogens, BT agents, moulds
- Easy access and easy GUI, result overview similar to the MBT software
- Comes with comprehensive microbiology knowledge base

MALDI Biotyper

Download available libraries
e.g. National Institute for Health

NIH provides a library to be downloaded
in addition to the Bruker MBT Libraries

Journal of
Clinical Microbiology

JCM Article | Journal Info. | Authors | Reviewers | Permissions | Journals.ASM.org

J Clin Microbiol. 2013 Mar; 51(3): 828–834. PMCID: PMC3592033
doi: [10.1128/JCM.02852-12](https://doi.org/10.1128/JCM.02852-12)

Development of a Clinically Comprehensive Database and a Simple Procedure for Identification of Molds from Solid Media by Matrix-Assisted Laser Desorption Ionization–Time of Flight Mass Spectrometry

[Anna F. Lau,](#)^a [Steven K. Drake,^b \[Leslie B. Calhoun,\]\(#\)^a \[Christina M. Henderson,\]\(#\)^a and \[Adrian M. Zelazny\]\(#\)^a](#)

[Author information](#) ► [Article notes](#) ► [Copyright and License information](#) ►

The MBT “**Open Microbiology**” concept enables this initiative!

MBT Subtyping Module

Listeria monocytogenes


- **Characteristic peaks** in the *Listeria* spectra are used to confirm the identification, using direct transfer sample preparation:

Without the MBT Subtyping Module the optimal way to achieve 100% ID is through extraction

With the MBT Subtyping Module 100% correct ID is obtained after direct transfer

Method	Accuracy for genus ID	Accuracy for species ID
Direct Transfer	100%	96.5%
Extended Direct Transfer	100%	95.4%
Formic Acid Extraction	100%	100%
Direct Transfer + MBT Subtyping Module	100%	100%

MBT Subtyping Module

Listeria identification


Result in MBT Compass

ID performed on direct smear => only if ID > 2.0 => further typing performed

NO additional steps required

ID	Name	Position	Chip	Detected Species	Score	Subtype
BTS	Bts	H12	0	Escherichia coli	2.45	
18141701	Sample	D1	0	Listeria innocua	2.18	typed as L. innocua
18141702	Sample	D2	0	Listeria innocua	2.02	typed as L. innocua
18141703	Sample	D3	0	Listeria innocua	2.21	typed as L. innocua
18141704	Sample	E1	0	Listeria monocytogenes	2.30	typed as L. monocytogenes
18141705	Sample	E2	0	Listeria monocytogenes	2.31	typed as L. monocytogenes
18141706	Sample	E3	0	Listeria monocytogenes	2.09	typed as L. monocytogenes
18141707	Sample	F3	0	Listeria welshimeri	2.18	typed as L. welshimeri
18141708	Sample	F4	0	Listeria welshimeri	2.21	typed as L. welshimeri
18141709	Sample	F5	0	Listeria welshimeri	2.05	typed as L. welshimeri
181417010	Sample	A1	0	Listeria innocua	2.11	presumptive L. monocytogenes, repeat measurement...
181417011	Sample	A2	0	Listeria ivanovii	2.25	typed as L. ivanovii
181417012	Sample	A3	0	Listeria ivanovii	2.23	typed as L. ivanovii
181417013	Sample	C1	0	Listeria seeligeri	2.10	calibrant/baseline too low, repeat measurement
181417014	Sample	C2	0	Listeria seeligeri	1.97	
181417015	Sample	C3	0	Listeria seeligeri	2.11	XXX typed as L. seeligeri


MBT Subtyping Module

Detection of specific antibiotic resistance

e.g. *S. aureus* typing for MRSA detection

Meticillin Resistens Staphylococcus Aureus

- The *S. aureus* subtyping application looks to detect the PSM-mec peptide peak present in a part of MRSA strains
 - If PSM-mec present: MRSA
(Very high % of correct detections)
 - If not, does not mean it is not a MRSA
- The proportion of MRSA with the PSM-mec peak is variable, depending on local epidemiology
 - Success rate varies from <10-90% depending upon region
- The *S. aureus* subtyping allows to detect immediately a part of MRSA strains, with no extra effort!

**Early Warning System
for Cattle & Raw Milk**

MALDI Biotyper

Official Method of Analysis

by AOAC International

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AOAC First Action *Official Method*SM 2017.09

Confirmation and Identification of
Salmonella spp.,
Cronobacter spp.


and other **gram-negative organisms**
by the Bruker MALDI Biotyper Method

AOAC First Action *Official Method*SM 2017.10

Confirmation and Identification of
Listeria monocytogenes,
Listeria spp.

and other **gram-positive organisms**
by the Bruker MALDI Biotyper method

AOAC-OMA: Association of Official Analytical
Chemists - Official Method of Analysis

MALDI Biotyper

ISO 16140-Part 6 Certification

by MicroVal

EU tanúsítvány

Certificate n° 2017LR73
Confirmation of *Salmonella* spp.
by the Bruker MALDI Biotyper method

Certificate n° 2017LR72
Confirmation of *Cronobacter* spp.
by the Bruker MALDI Biotyper method

Certificate n° 2017LR75
Confirmation of *Listeria* spp. and *Listeria monocytogenes*
by the Bruker MALDI Biotyper method

CERTIFICATE OF COMPLIANCE

LLOYD'S REGISTER QUALITY ASSURANCE

hereby declares that the certification assessment has demonstrated that

MALDI Biotyper®

Complete Solution for the confirmation of *Listeria* spp. and *Listeria monocytogenes*

Manufactured and supplied by:
Bruker Daltonik GmbH
Fahrenheitstraße 4
D-28359 Bremen
GERMANY

has been validated and revealed to be at least equivalent to the reference method as demonstrated by the validation study report. The summary of the validation report is available on the MicroVal website: www.microval.org

Reference methods:

1. ISO 11290-1 (2017): Microbiology of the food chain - Horizontal method for the detection and enumeration of *Listeria monocytogenes* and other *Listeria* spp. - Part 1: detection method;
2. ISO 11290-2 (2017): Microbiology of the food chain - Horizontal method for the detection and enumeration of *Listeria monocytogenes* and other *Listeria* spp. - Part 2: enumeration method.

Scope: Confirmation of *Listeria* spp. and *Listeria monocytogenes* from colonies isolated on Ottaviani & Agosti Agar, Palcam, Oxford, Modified Oxford and RAPID'L mono, and any non-selective nutrient agars

The validation and certification has been performed in accordance with ISO/DIS 16140-6:2017 and the MicroVal Rules and Certification Scheme version 8.

Certificate no.: 2017LR75

First approval date: 26 April 2018
Expiry date: 25 April 2022

ISSUED BY: Lloyd's Register Nederland B.V.,
Rotterdam, The Netherlands

Certificate no.: 2017LR75

03-05-2018

Page 1 of 1

K.P. v.d. Mandelsteilaan 41a, 3062 MB Rotterdam, The Netherlands. KvK nr. 24247948
This approval is carried out in accordance with the LRQA assessment and certification procedures and monitored by LRQA

Lloyd's Register
an Lloyd's
list company

Lloyd's Register
an Lloyd's
list company

Lloyd's Register Group Limited, its affiliates and subsidiaries, including Lloyd's Register Quality Assurance Limited (LRQA), and their respective officers, employees or agents are, individually and collectively, referred to in this document as 'Lloyd's Register'. Lloyd's Register assumes no responsibility and shall not be liable to any person for any loss, damage or expense caused by reliance on the information or advice in this document or for any other reason, unless that person has signed a contract with the relevant Lloyd's Register entity for the provision of this information or advice and in that case any responsibility or liability is excluded to the extent permitted by law.

Confirmation workflow: comparison

Biochemical tests – MALDI Biotyper

Methods	Analysis step	<i>Listeria</i>	<i>L. monocytogenes</i>	<i>Salmonella</i>	<i>Cronobacter</i>
Reference FDA, USDA, ISO methods	Purification step on a non-selective agar (+ 24 hrs)	• Required	• Required	• Required except in the ISO 6579	• Required
	Confirmation of characteristic colonies (up to 5)	<ul style="list-style-type: none"> • Gram staining • Catalase • (Motility) 	<ul style="list-style-type: none"> • Gram staining • Catalase • Hemolysis • Biochemical tests • (Motility) 	<ul style="list-style-type: none"> • Biochemical tests • Serological tests 	<ul style="list-style-type: none"> • Oxidase tests • Biochemical tests
	Final result if presence of characteristic colonies on selective media	Day 1	Day 2	Day 2 to 3	Day 2
	Handling time*	Not available	> 10 min	4 min	3,5 min
MBT	Direct Transfer from isolated colonies on selective agars	<ul style="list-style-type: none"> • No purification step required • Same workflow for all the pathogens • Result in 15 min to 1 hour (1 to 95 isolates) • 20 sec handling time / isolate 			

*handling time from ISO 16140 Validation reports


Workflow: comparison

16S rDNA sequencing – MALDI Biotyper


Mikrobiológiai teszt

16S rDNA workflow

MBT workflow

MALDI Biotyper

Basic Bruker consumables

Quality controlled

- MALDI matrix:
HCCA Matrix
- Standard for QC:
Bacterial Test Standard
- Disposable MALDI target:
MBT Biotarget 96

Compared to traditional ID methods,
significantly reduced consumables budget
needed for standard ID!

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MALDI Biotyper

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ID of Gram +/- bacteria and yeast (2748 species covered)


Special library for filamentous fungi (152 species/species groups)

Automated calibration and QC by quality controlled BTS

Completely usable disposable MALDI targets

Paperless & traceable workflow

Sample prep facilitated by MBT Pilot and MBT Galaxy

www.bruker.com